



Inactivated *Lactobacillus reuteri* Culture Fluid (DAG-WT3693)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	Control
Bio-activity	$\geq 10^6$ CFU/mL
Format	Liquid
Size	1 mL
Buffer	DMEM
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	<i>Lactobacillus reuteri</i> ; <i>L. reuteri</i>
----------	--