



Inactivated *Cryptococcus neoformans* Culture Fluid (DAG-WT3697)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	Control
Bio-activity	≥10e6 CFU/mL
Format	Liquid
Size	1 mL
Buffer	DMEM
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	<i>Cryptococcus neoformans</i> ; <i>C. neoformans</i>
----------	---