



# Pentadecafluorooctanoic acid [BSA] (DAG-WT3707B)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

## BACKGROUND

Keywords	335-67-1; Pentadecafluorooctanoic acid; PFOA; Perfluorocaprylic acid
----------	--