



Tetrachlorobiphenyl [OVA] (DAG-WT3708O)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	2,3,4,5-Tetrachlorobiphenyl; 33284-53-6; 1,1'-Biphenyl, 2,3,4,5-tetrachloro-; PCB 61
----------	--