



Amitraz [BSA] (DAG-WT3711B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Amitraz; 33089-61-1; 1,5-Di(2,4-dimethylphenyl)-3-methyl-1,3,5-triazapenta-1,4-diene
----------	--