



Cyprodinil [BSA] (DAG-WT3722B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Cyprodinil; 121552-61-2; 4-Cyclopropyl-6-methyl-N-phenylpyrimidin-2-amine; 4-Cyclopropyl-6-methyl-N-phenyl-2-pyrimidinamine
----------	---