



Isofetamid [BSA] (DAG-WT3724B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Isofetamid; 875915-78-9; IKF 5411; IKF-5411; UNII-95499Z1H0H
----------	--