



Tebufenozide [BSA] (DAG-WT3728B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords Tebufenozide; 112410-23-8; N-tert-Butyl-N'-(4-ethylbenzoyl)-3,5-dimethylbenzohydrazide