



Quinoxifen [OVA] (DAG-WT3730O)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Quinoxifen; 124495-18-7; 5,7-dichloro-4-(4-fluorophenoxy)quinoline; Quinoline, 5,7-dichloro-4-(4-fluorophenoxy)-; DE 795
----------	--