



Bupirimate [BSA] (DAG-WT3732B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|--------------|------------------|
| Conjugate | BSA |
| Applications | ELISA, LFIA |
| Format | Liquid |
| Size | 1 mg |
| Buffer | 10mM PBS, pH 7.4 |
| Preservative | None |
| Storage | Store at -20°C |
| Ship | Wet ice |

BACKGROUND

| | |
|----------|---|
| Keywords | Bupirimate; 41483-43-6; Nimrod; 5-Butyl-2-ethylamino-6-methylpyrimidin-4-yl dimethylsulfamate |
|----------|---|