



Diethofencarb [BSA] (DAG-WT3737B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Diethofencarb; 87130-20-9; powmil; Isopropyl 3,4-diethoxycarbanilate
----------	--