



Thiachloprid [BSA] (DAG-WT3740B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Thiachloprid; 111988-49-9; Calypso; (3-((6-Chloro-3-pyridinyl)methyl)-2-thiazolidinylidene)cyanamide; (Z)-thiachloprid
----------	--