



Fenazaquin [BSA] (DAG-WT3741B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Fenazaquin; 120928-09-8; Magister; 4-(4-(tert-Butyl)phenethoxy)quinazoline
----------	--