



Nitroxinil [BSA] (DAG-WT3742B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Nitroxynil; Nitroxinil; 1689-89-0; 4-Hydroxy-3-iodo-5-nitrobenzonitrile; Dovenix
----------	--