



Sulfometuron [BSA] (DAG-WT3748B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Sulfometuron; 74223-56-6; 2-[(4,6-dimethylpyrimidin-2-yl)carbamoylsulfamoyl]benzoic acid; 2-(3-(4,6-Dimethylpyrimidin-2-yl)ureidosulfonyl)benzoic acid
----------	--