



Thidiazuron [BSA] (DAG-WT3749B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Thidiazuron; 51707-55-2; 1-Phenyl-3-(1,2,3-thiadiazol-5-yl)urea Urea, N-phenyl-N'-1,2,3-thiadiazol-5-yl-; 1-phenyl-3-(thiadiazol-5-yl)urea
----------	---