



## Cinosulfuron [OVA] (DAG-WT3750O)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|              |                  |
|--------------|------------------|
| Conjugate    | OVA              |
| Applications | ELISA, LFIA      |
| Format       | Liquid           |
| Size         | 1 mg             |
| Buffer       | 10mM PBS, pH 7.4 |
| Preservative | None             |
| Storage      | Store at -20°C   |
| Ship         | Wet ice          |

### BACKGROUND

|          |                                                                                                                      |
|----------|----------------------------------------------------------------------------------------------------------------------|
| Keywords | Cinosulfuron; 94593-91-6; CGA 142464; 1-(4,6-dimethoxy-1,3,5-triazin-2-yl)-3-[2-(2-methoxyethoxy)phenyl]sulfonylurea |
|----------|----------------------------------------------------------------------------------------------------------------------|