



Bensulfuron-methyl [OVA] (DAG-WT3751O)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Bensulfuron-methyl; Bensulfuron methyl; 83055-99-6; Londax; Methyl bensulfuron
----------	--