



Fenoxaprop [OVA] (DAG-WT3757O)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Fenoxaprop; 95617-09-7; Fenoxaprop acid; 2-[4-[(6-chloro-1,3-benzoxazol-2-yl)oxy]phenoxy]propanoic acid
----------	---