



Dicamba [BSA] (DAG-WT3759B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Dicamba; 1918-00-9; 3,6-Dichloro-2-methoxybenzoic acid; Benzoic acid, 3,6-dichloro-2-methoxy-; Banvel
----------	---