



Bifenox [BSA] (DAG-WT3761B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Bifenox; 42576-02-3; Methyl 5-(2,4-dichlorophenoxy)-2-nitrobenzoate; Modown; 5-(2,4-Dichlorophenoxy)-2-nitrobenzoic acid methyl ester
----------	---