



Picloram [BSA] (DAG-WT3766B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Picloram; Tordon; Grazon; 4-Amino-3,5,6-trichloropyridine-2-carboxylic acid; 1918-02-1
----------	--