



Chloridazon [BSA] (DAG-WT3768B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Chloridazon; 5-Amino-4-chloro-2-phenylpyridazin-3(2H)-one; Chloridazone; 1698-60-8
----------	--