



Recombinant Hepatitis E Virus (HEV) VLP (DAG-WT377)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	25 µl, 50 µl
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -80°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Introduction	Hepatitis E is a liver infection caused by the hepatitis E virus (HEV). HEV is found in the stool of an infected person. It is spread when someone unknowingly ingests the virus – even in microscopic amounts.
Keywords	Hepatitis E Virus; HEV; VLP