



2-hydroxy-3,5-dinitrobenzohydrazide [BSA] (DAG-WT3770B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords 955-07-7; 2-hydroxy-3,5-dinitrobenzohydrazide; 3,5-dinitrosalicylhydrazide; Benzoic acid, 2-hydroxy-3,5-dinitro-, hydrazide; 3,5-Dinitrosalicylohydrazide