



Nimesulide [BSA] (DAG-WT3773B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	nimesulide; 51803-78-2; N-(4-Nitro-2-phenoxyphenyl)methanesulfonamide; Mesulid
----------	--