



## Fenamic acid [BSA] (DAG-WT3782B)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|              |                  |
|--------------|------------------|
| Conjugate    | BSA              |
| Applications | ELISA, LFIA      |
| Format       | Liquid           |
| Size         | 1 mg             |
| Buffer       | 10mM PBS, pH 7.4 |
| Preservative | None             |
| Storage      | Store at -20°C   |
| Ship         | Wet ice          |

### BACKGROUND

|          |                                                                              |
|----------|------------------------------------------------------------------------------|
| Keywords | N-phenylanthranilic acid; 91-40-7; Fenamic acid; 2-(Phenylamino)benzoic acid |
|----------|------------------------------------------------------------------------------|