



Fenamic acid [OVA] (DAG-WT3782O)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	N-phenylanthranilic acid; 91-40-7; Fenamic acid; 2-(Phenylamino)benzoic acid
----------	--