



Bupivacaine [BSA] (DAG-WT3784B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords bupivacaine; 2180-92-9; 1-butyl-N-(2,6-dimethylphenyl)piperidine-2-carboxamide; 38396-39-3