



## Benzafibrate [BSA] (DAG-WT3786B)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

### BACKGROUND

Keywords	Bezafibrate; 41859-67-0; Bezalip; Cedur; Bezafibrat
----------	---