



Bisoxatin [BSA] (DAG-WT3790B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	Bisoxatin; 17692-24-9; 2,2-bis(4-hydroxyphenyl)-4H-1,4-benzoxazin-3-one; 2,2-Bis(p-hydroxyphenyl)-2H-1,4-benzoxazin-3(4H)-one
----------	---