



Butyl hydroxyanisole [BSA] (DAG-WT3792B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords Butyl hydroxyanisole; 2-tert-Butyl-4-methoxyphenol; 3-tert-Butyl-4-hydroxyanisole; 121-00-6; 25013-16-5; 4-Hydroxy-3-tert-butylanisole