



2,4-toluylene diisocyanate [OVA] (DAG-WT3797O)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Keywords	584-84-9; 2,4-Diisocyanato-1-methylbenzene; 2,4-Toluene diisocyanate; Toluene-2,4-diisocyanate
----------	--