



# Recombinant *Campylobacter jejuni* CdtB protein [His] (DAG-WT3805)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA, WB
Molecular Weight	27.5 kDa
Format	Liquid
Size	1 mg
Buffer	20mM Tris-HCl, pH 8.0, 0.2M NaCl, 2mM DTT
Preservative	None
Storage	Store at -80°C. Avoid multiple freeze/thaw cycles.
Ship	Dry ice

## BACKGROUND

Keywords	<i>Campylobacter jejuni</i> ; C. Jejuni; CdtB; Cytolethal distending toxin B
----------	--