



# Recombinant Hantaan Virus Nucleocapsid (N) Protein [His] (DAG-WT3806)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Purity</b>	> 95% as determined by SDS-PAGE
<b>Conjugate</b>	His
<b>Applications</b>	ELISA
<b>Format</b>	Liquid
<b>Size</b>	100 µg, 500 µg, 1 mg
<b>Buffer</b>	PBS, pH 7.4
<b>Preservative</b>	0.05% Sodium azide
<b>Storage</b>	Store at -20°C

## BACKGROUND

<b>Keywords</b>	Hantaan Virus; Nucleocapsid; HTNV
-----------------	-----------------------------------