



# Recombinant Human Thrombomodulin (DAG-WT3811)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	Calibration, Standard
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -20°C

## BACKGROUND

Keywords	Thrombomodulin; TM; CD141; BDCA-3
----------	-----------------------------------