



Recombinant Human P-selectin protein (DAG-WT3812)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
---------------	---------------------------------

Conjugate	His
------------------	-----

Applications	Calibration, Standard
---------------------	-----------------------

Format	Liquid
---------------	--------

Size	1 mg
-------------	------

Buffer	PBS, pH 7.4
---------------	-------------

Preservative	None
---------------------	------

Storage	Store at -20°C
----------------	----------------

BACKGROUND

Keywords	P-selectin; SELP; Cell adhesion molecule; CAM
-----------------	---