



# Recombinant Human AAG protein (DAG-WT3813)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>Conjugate</b>	His
<b>Applications</b>	Calibration, Standard
<b>Format</b>	Liquid
<b>Size</b>	1 mg
<b>Buffer</b>	PBS, pH 7.4
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C

## BACKGROUND

<b>Keywords</b>	AAG; alpha-1-acid glycoprotein; $\alpha$ 1AGp; AGP
-----------------	--