



Recombinant Canine Distemper Virus F Antigen (DAG-WT3817)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Canine
Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	Calibration, Standard
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	CDV; Canine Distemper virus
----------	-----------------------------