



# Recombinant Feline pancreatic lipase (PL) (DAG-WT3825)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|              |                                 |
|--------------|---------------------------------|
| Species      | Feline                          |
| Purity       | > 90% as determined by SDS-PAGE |
| Conjugate    | His                             |
| Applications | Calibration, Standard           |
| Format       | Liquid                          |
| Size         | 1 mg                            |
| Buffer       | PBS, pH 7.4                     |
| Preservative | None                            |
| Storage      | Store at -20°C                  |

## BACKGROUND

|          |                       |
|----------|-----------------------|
| Keywords | Pancreatic lipase; PL |
|----------|-----------------------|