



Recombinant Human MS4A6A protein [GST] (DAG-WT3829)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Conjugate	GST
Applications	ELISA, WB
Molecular Weight	52.8 kDa
Format	Liquid
Size	10 µg, 25 µg
Buffer	50 mM Tris-HCl, 10 mM reduced Glutathione, pH8.0
Preservative	None
Storage	Store at -80°C

BACKGROUND

Keywords	MS4A6A; Membrane spanning 4-domains A6A
----------	---

GENE INFORMATION

UniProt ID	AAH22854
------------	----------