



Recombinant Human TRIM33 protein (DAG-WT3830)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Conjugate	TBD
Applications	ELISA, WB
Molecular Weight	124 kDa
Format	Liquid
Size	1 mg
Buffer	50 mM Tris-HCl, 10 mM reduced Glutathione, pH8.0
Preservative	None
Storage	Store at -80°C

BACKGROUND

Keywords	TRIM33; TIF1 gamma; TIF1-γ
----------	----------------------------

GENE INFORMATION

UniProt ID	Q9UPN9
------------	--------