



Recombinant CMV p38-pp65-gp52-pp150 chimeric protein (DAG-WT3833)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Format	Liquid
Size	1 mg
Buffer	20mM Tris-HCl, pH 8.0, 100mM NaCl, 0.5mM EDTA, 0.15mM DTT
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	Cytomegalovirus; CMV; HCMV; HHV-5; p38; pp65; gp52; pp150
----------	---