



Recombinant Kunjin virus NS1 [His] (DAG-WT3838)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	WB, ELISA
Format	Liquid
Size	1 mg
Buffer	Tris/PBS-based buffer
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	Kunjin virus; KUNV; Non-structural protein 1; NS1
----------	---