



# Recombinant Myxovirus resistance protein 1 (MxA) [His] (DAG-WT3839)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Species</b>	Human
<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>Conjugate</b>	His
<b>Applications</b>	WB, ELISA
<b>Molecular Weight</b>	76.3 kDa
<b>Reconstitution</b>	Add deionized water to prepare a working stock solution of approximately 0.5 mg/ml and let the lyophilized pellet dissolve completely at 37°C
<b>Format</b>	Lyophilized
<b>Size</b>	100 µg
<b>Buffer</b>	20 mM Tris buffer, 50 mM NaCl, 5% (w/v) trehalose, pH 7.5
<b>Preservative</b>	None
<b>Storage</b>	Store at -80°C

## BACKGROUND

**Keywords** Interferon-induced GTP-binding protein Mx1; MX1; MxA; Myxovirus resistance protein 1