



Recombinant Middle East Respiratory Syndrome-CoV (MERS-CoV) VLP (DAG-WT384)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -80°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Introduction	MERS-CoV is a zoonotic virus, which means it is a virus that is transmitted between animals and people. Studies have shown that humans are infected through direct or indirect contact with infected dromedary camels. MERS-CoV has been identified in dromedaries in several countries in the Middle East, Africa and South Asia.
Keywords	Middle east respiratory syndrome-CoV; MERS-CoV; VLP