



## NP [AECM-FICOLL] (DAG-WT3840)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	AECM-FICOLL
Applications	ELISA
Format	Liquid
Size	10 mg, 100 mg
Buffer	PBS
Preservative	None
Storage	Store at -20°C

### BACKGROUND

Keywords	NP; 4-Hydroxy-3-nitrophenylacetic; AminoEthylCarboxyMethyl; FICOLL
----------	--