



Recombinant *Brucella abortus* Omp28 [His] (DAG-WT3844)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA, WB
Format	Liquid
Size	1 mg
Buffer	Tris-based buffer
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	<i>Brucella abortus</i> ; Outer membrane protein; Omp28
----------	---