



Native Human Neutrophil Elastase (DAG-WT3845)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	N/A
Reconstitution	Reconstitute in 50 mM sodium acetate, pH 5.5, with 200 mM NaCl
Bio-activity	≥20 units/mg
Format	Lyophilized
Size	10 µg, 100 µg
Storage	Store at -20°C

BACKGROUND

Keywords	Elastase; EC 3.4. 21.36; Serine protease
----------	--