



Recombinant Human Hyaluronidase (DAG-WT3846)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	47 kDa
Format	Liquid
Size	10 µg, 1 mg
Buffer	PBS with 10% glycerol
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	Hyaluronidase; 3.2.1.35; HYAL
----------	-------------------------------